step-by-step **Easy Orthotic Solutions** tutorials for you to try for the Fingers, Hand, **Thumb, Wrist and Forearm**

A booklet of easy-to-follow orthotic instructions using Orficast and Orficast More



IN ORTHOTIC FABRICATION FOR PHYSICAL REHABILITATION

today!



Orficast is a unique thermoplastic material based on threads of thermoplastic. It is supplied on a roll which makes it easy to size, cut and store as opposed to most splinting materials.

Why therapists **love Orficast**

- in your workplace.
- with regular scissors.
- ventilation of the skin.

→ The compact product size takes up little space

 \rightarrow A convenient dispenser box lets you cut only the required amount of material that you need. The rest of the product will remain neat and tidy in the box.

→ Orficast and Orficast More are easy to cut

 \rightarrow The airy, knitted structure allows for

→ This product increases comfort on patients with fragile and/or sensitive skin.

 \rightarrow Orficast molds easily to the shape of the body part allowing for a snug fit to prevent pressure points.

→ When heated, Orficast and Orficast More stick firmly to itself, but not to skin or hair.

The Product Line

Orficast

Orficast in 3 cm and 6 cm widths can be used for all types of finger and hand-based immobilization orthoses.

Length	3 meters (around 10 feet)
Widths	3 cm (1") and 6 cm (2")
Colors	Blue, Black and Orange



Orficast More

Orficast More in 6 cm, 15 cm and 30 cm widths is a thicker product that offers increased rigidity, and more support. Orficast More can be used for larger orthoses, such as thumb, wrist and forearm-based orthoses.

Length	3 meters (around 10 feet)
Widths	6 cm (2"), 15 cm (6") and 30 cm (12")
Colors	Blue, Black and Orange



Questions?

We invite you to explore Orficast and Orficast More's versatility, ease of use and working properties by going through the step-by-step instructions in this booklet.

Do you have a question about Orficast?

Do you need help with an orthotic fabrication challenge?

Email your question to welcome@orfit.com.

Our product specialists will be happy to help you.



Table of Contents

This booklet shares ideas and fabrication techniques for a variety of the most common orthoses made with Orficast and Orficast More.

Join with thousands of therapists who utilize Orficast and Orficast More for orthotic fabrication.

Area of Focus

Orthotic Solution

Mallet finger	Mallet Orthosis Ribbon De
Mallet finger	Mallet Orthosis Crossover
Boutonniere deformity	Finger Wrap Extension Or
Proximal Interphalangeal (PIP) joint injury	Buddy Tapes
Swan Neck Deformity	Anti-Swan Neck Orthosis
Trigger finger	Trigger Finger Orthosis
Metacarpal Phalangeal (MCP) joint sprain	Relative Motion Orthosis
Claw hand from ulnar nerve palsy	Anti-Claw Orthosis
Thumb sprain, pain, arthritis or similar	Short Thumb Opponens O
Wrist sprain, tendinitis, fracture or similar	Wrist Orthosis - Thumb H
Wrist and thumb pain, arthritis, sprain, tendinitis or similar	Long Opponens Orthosis
Forearm bones, distal radial-ulnar joint and/or TFCC ligament	Muenster Orthosis

Material

Page

esign	Orficast 3 cm (1")	14
er Design	Orficast 3 cm (1")	16
rthosis	Orficast 3 cm (1")	18
	Orficast 3 cm (1")	20
	Orficast 3 cm (1")	22
	Orficast 3 cm (1")	24
	Orficast 6 cm (2")	26
	Orficast More 6 cm (2")	28
Orthosis	Orficast More 15 cm (6")	30
lole Design	Orficast More 15 cm (6") and Orficast More 6 cm (2")	32
	Orficast More 15 cm (6")	34
	Orficast More 30 cm (12") and Orficast More 6 cm (2")	36

browse this booklet to start your orthotic fabrication





Tips and Tricks

for Working with Orficast and Orficast More

How to use Orficast to create How to clean an Orficast finger cuffs with thread. Orthosis.

• Use Orficast 3 cm (1") to create a finger cuff for a mobilization orthosis. Simply activate a small piece of Orficast and while it is warm, fold it over the elastic thread or static line and mold around your patient's finger.

• If you only have the widest version of Orficast More in 30 cm (12") you can easily cut it lengthwise along the thermoplastic threads to make Orficast More in other widths. It will not unravel.

Orficast and Orficast More orthoses can be easily cleaned in the washing machine if placed in a lingerie bag and using a cold-water cycle. If the straps are adhered well to the orthosis, they will not come off.

Abbreviations

MCP	Metacarpal-pha
PIP	Proximal interpl
DIP	Distal interphala

How to use Orficast to fabricate ring style orthoses like the anti-swan neck orthosis?

- Orficast can be rolled into a thin tubular shape or folded in half to fabricate ring style orthoses.
- Use Orficast rolled into tubular shapes to create finger separators in resting hand orthoses.
- Simply use dry heat from a heat gun on the volar portion of the resting hand orthosis and on the Orficast to adhere together.
- Always check the fit of the orthosis after fabrication to avoid any pressure points.



See our Orfit Academy courses on Finger Orthoses for more fabrication instructions.

alangeal joint

- phalangeal joint
- langeal joint



How to attach straps with the help of Edging Strip

Here is a new, quick and easy method to attach hook and loop (e.g. Velcro) straps to Orficast. We will be using an unexpected ally for an unbreakable bond: Thermoplastic Edging Strip.

The Edging Strip is a roll of super-thin low-temperature thermoplastic material, traditionally used to line the edges of an orthosis. This strip material will bond quickly to Orficast, even if it's still damp. A perfect tool to attach straps to Orficast!

You will need:

- Orfit edging strip
- Loop strap
- Adhesive-backed hook
- Orficast scraps
- Heat gun
- Scissors











1. Prepare a piece of adhesive-backed hook and a similar-sized piece of Orfit Edging Strip.



2. Heat up the adhesive glue side of the hook briefly with the heat gun and attach to the Orfit Edging Strip. Trim the corners.







3. Place this combination in the water bath and let the Edging Strip become transparent. Once ready, attach the hook strap firmly to the molded orthosis.



4. Briefly heat up the other side of the orthosis with the heat gun and press the fuzzy side of the loop strap directly into the heated thermoplastic.







5. Use the heat gun to heat up a small scrap of Orficast material and place it on top of the loop strap to reinforce its attachment. If necessary, trim the excess scrap material for a neat finish.



Materials	Preparation	
 Orficast 3 cm (1"). 2.5 cm (0.5") adhesive backed hook and loop strapping. 	1. Cut a strip about 15 (6-8") long.	
	2. Activate by immers 65°C (160°F).	
	 Pat dry on towel br material in half so t 7.5-10 cm (3-4") lon 	

Fabricate



- 1. Place middle portion underneath the fingertip.
- 2. Bring one side of the material to meet the PIP joint dorsally.



- 5. Let the material harden with the fingertip held in extension or slight hyperextension.
- 4. Smooth the layers together.

3. Wrap the other

side of the material

overlapping on the

around the finger

to the PIP joint

first layer.


5-20 cm

sion in water bath

riefly and fold that piece is now

Position

- Maintain DIP joint extension.
- Pearl: Strap snugly over middle phalanx to keep orthosis in place.
- Check patient's ability to hold DIP extension against gravity.







6. Remove the orthosis, 7. Heat the dorsal trim the corners and the sides.



portion of the orthosis with the heat gun and adhere the loop strap and hook by pressing firmly into the Orficast material. Strap snugly around the middle phalanx.









Materials

- Orficast 3 cm (1").
- ✓ 2.5 cm (0.5") adhesive backed hook and loop strapping.

Preparation

- 1. Cut a strip about 15-20 cm (6-8") long.
- 2. Activate by immersion in water bath 65°C (160°F).
- 3. Pat dry on towel briefly and fold material in half so that piece is now 7.5-10 cm (3-4") long.

Fabricate



1. Fold Orficast in 2. Place on dorsum half and cut up the middle leaving about 2 cm (1") uncut.



of finger with uncut

edge near PIP joint.

3. Wrap one wing under 4. Bond well by finger tip on one side and the other wing under finger tip on the other side.

Position

- Maintain DIP joint extension
- Pearl: Strap snugly over middle phalanx to keep orthosis in place.
- Check patient's ability to hold DIP extension against gravity.



Wearing schedule is full time for 6-8 weeks.



pressing together as you hold the fingertip in extension until the material hardens.



5. Place a strap over the middle phalanx.



Finger Wrap Extension Orthosis

CLINICAL CONDITION Boutonniere Deformity

Materials

V Orficast 3 cm (1").

Preparation

- 1. Cut one strip about 20 cm (8") long.
- 2. Activate by immersion in water bath 65°C (160°F).
- 3. Pat dry with towel briefly.

Position

- Keep PIP joint in maximum extension.

- and cannot be pulled apart.

Fabricate



- 1. Wrap one end of the 2. Continue to wrap Orficast material around the distal phalanx of the finger and secure together.
- around the finger on a diagonal and with a slight stretch.

18

Avoid putting pressure on the dorsal PIP joint.


• Pearl: Complete the orthosis and then address positioning.

• Make sure all layers are smoothed together.

• Make sure that the distal edge of Orficast is adhered to the orthosis









4. Rub all the layers together to smooth them and position the finger as needed. The orthosis can be removed for trimming or to dry heat and secure a corner that did not stick.

Buddy Tapes

CLINICAL CONDITION Proximal Interphalangeal Joint Injury



V Orficast 3 cm (1").


Preparation

- 1. Cut two strips about 10 cm (5") long.
- 2. Activate by immersion in a water bath 65°C (160°F).
- 3. Pat dry with a towel and either roll into a tube shape or fold in half lengthwise.
- 4. Maintain the fist in MCP flexion which makes sure that the PIP joints are at the same level while molding.

Fabricate



1. Wrap the material over the palmar surface of the proximal phalanx and connect it to itself as a ring with pressure.



Position

3. Repeat the same for 2. Continue wrapping the second part of around the proximal the orthosis which second finger and rests over the middle pass it around to phalanges. connect on the volar surface. Connect by

phalanx of the

applying sufficient pressure. Trim excess away.

• Buddy tapes are placed over proximal and middle phalanges

• Pearl: Do not block PIP or DIP joint motion

• Position the MCP joints in 45 degrees of flexion and PIP and DIP joints in extension so that the proximal phalanges are aligned.







. Make sure the buddy tapes can be removed 4 easily. Dip in hot water and stretch if needed to enlarge each ring.



V Orficast 3 cm (1").

Preparation

- 1. Cut one strip about 20 cm (8") long.
- 2. Activate by immersion in a water bath 65°C (160°F).
- 3. Pat dry with a towel and either roll into a tube shape or fold in half lengthwise.
- 4. Maintain the fist in MCP flexion which makes sure that the PIP joints are at the same level while molding.

Fabricate



1. Roll the Orficast into a tight tube.



- 2. Mold an oval shape over the dorsum of the first phalanx and the middle phalanx, keeping the finger in PIP joint flexion. snugly.
- 3. Begin at the PIP joint on one side and complete the oval in the same

Position

• Orthosis will block 10-15 degrees of end range PIP joint extension.

• Pearl: Keep PIP joint in maximum flexion during fabrication.



spot, attaching the material to itself



4. Slip the remaining material underneath the flexed PIP joint and attach to the opposite side of the oval. Pinch the material together underneath the PIP joint.



5. Trim the remaining material away and hold in PIP joint flexion until hardened.



- 1. Cut one strip about 20 cm (8") long.
- 2. Activate by immersion in water bath 65°C
- 3. Pat dry with towel briefly.

Position



1. Wrap one end of the 2. Fold the material on Orficast material around the proximal phalanx of the finger. Seal it together.



a diagonal so that

the long end extends

the material from the

palm and fold it back

on itself towards the

finger all the way up to the PIP joint

creating a double

layer.

into the palm. Take

3. Now fold over again from the PIP joint back to the palm to create another layer of material on top. There is now one layer around the proximal phalanx, and three layers on top of each other in front of the MCP crease into the palm.

• Keep MCP joint hyperextended during fabrication

• Pearl: Make sure finger does not trigger with orthosis. If yes, block flexion of the PIP joint as well.





4. Let harden, then remove and trim the corners.

Relative Motion Extension Orthosis

CLINICAL CONDITION MCP Joint Sprain

Materials

Orficast 6 cm (2").

Preparation

- 1. Cut one strip about 20 cm (8") long or slightly wider than the patient's hand.
- 2. Activate by immersion in water bath 65°C (160°F).
- 3. Pat dry on towel briefly and fold

Position

Fabricate



1. Measure a piece of Orficast slightly wider than all four fingers.



2. Fold into thirds.



3. Place underneath the proximal phalanx of the middle finger.



4. Wrap one side around the index finger at the level of the proximal phalanx. Wrap the other side around the ring finger.

5. With the patient's hand in neighboring fingers.

• Make sure the intended MCP joint is held in more extension or more flexion than adjacent MCP joints.

• Pearl: Make sure the orthosis can be easily put on and off the fingers.

• Mark the orthosis so the patient knows exactly how to put it on.



supination, push into the middle finger with your thumb over the material. Keep the ring and index fingers relatively close to each other. Apply downward pressure over the middle finger so that the finger is more extended than its



6. Let harden and pronate the hand to observe. The middle finger should sit in slightly more extension at the MCP joint than the ring and index fingers. If needed, dip each ring into the splint pan and stretch so that the orthosis can be easily placed on and off the hand. Mark an arrow where the injured finger needs to rest.



• Make sure the MCP joints are in 60 degrees of flexion.

• Roll the material on the volar surface into a tight tubular shape

• Indent slightly on the dorsal surface to keep the fingers separated.

• Mark the orthosis so the patient knows exactly how to put it on.





4. Trim the excess material. Cut open and apply loop and hook strapping or leave closed for a slip-on orthosis.



5. Indent slightly between the fingers.

Short Thumb Opponens Orthosis

CLINICAL CONDITION Thumb pain, sprain, arthritis

Materials

- Orficast More 15 cm (6").
- Adhesive backed hook and loop strapping.

Preparation

1. Cut a 15 by 15 cm (6" by 6") square of material.

2. Cut off the corner at one end and into the center from the opposite corner, leaving enough width of the center portion to cover the web space between the IP crease and the distal palmar crease.

together at the

thumb IP crease and

all the way proximal

to the corners below


ulnar wings together


the wrist. Pinch the

with stretch.

3. Activate by immersion in water bath 65°C (160°F).

Fabricate





1. Place the center 2. Pinch the corners portion of the material into the web space and stretch around the radial side of the thumb.



3. While warm, pull the material away from the thumb and cut, forming a flattened seam. Open the ulnar pinch, and re-pinch until hardened.



Position

- Put the thumb in a functional position.
- Pearl: Make sure the IP joint is free; trim all rough edges away and reinforce the straps with additional scraps of Orficast.





4. Remove the orthosis 5. Attach a single strap and trim the corners and sides.



over the ulnar side.



Wrist Orthosis **Thumb Hole Design**

CLINICAL CONDITION Wrist pain, sprain, tendinitis, fractures and/or arthritis

Fabricate



1. Cut a small hole for the thumb, slightly off center at one end of the material. Place the thumb through the hole and stretch the material around the palm, pinching together over the metacarpals.



Materials

Orficast More 15 cm (6").

Orficast More 6 cm (2").

2. Stretch the material along the volar forearm and pinch on the dorsal surface proximally.



material harden.

3. Pinch once more dorsally at wrist level. Place in thumb hole and desired wrist sides. position, and let the

- Preparation
- 1. Cut a rectangle from Orficast More 15 cm (6") about 30 cm (10") long and a similar length of the Orficast More 6 cm (2").
- 2. Activate by immersion in water bath 65°C (160°F).
- 3. Press both materials together on a towel while patting excess water out. Compress by rolling in the towel.

Position

- The wrist is typically held in 20-30 degrees of extension.
- Pearl: Clear the thumb hole for comfort and allow for full flexion of all fingers at the distal palmar crease.







- and trim the corners.
- 4. Remove the orthosis 5. Apply straps.





Long Opponens Orthosis

CLINICAL CONDITION Thumb and wrist pain, sprains, arthritis and/or similar



Fabricate




1. Measure a rectangle of Orficast More 15 cm (6") long enough to cover the thumb and 2/3rd of the forearm.



2. Cut off the two corners, activate the material in hot water and pat dry.



3. Fold over the distal edge and place over the radial border of the thumb and wrist.



• The thumb is held in a functional position and the wrist is in 20-30 degrees of extension.

• Pearl: If the thumb IP joint is enlarged, make sure the overlapped portion around the thumb portion allows for easy on and off. If not consider a loop strapping closure instead of overlapping







4. Overlap the material 5. Attach straps.



See page 12: How to attach straps to wet Orficast and Orficast More.

35



Muenster Orthosis

CLINICAL CONDITION Prevention or Limitation of **Forearm Motion**

Materials

- Orficast More 6 cm (2").
- ✓ Orficast More 30 cm (12")
- . Adhesive backed hook and loop strapping.

Preparation

- 1. Cut a rectangle of the wider material about ³/₄ the length of the forearm and a strip of the narrow Orficast More the same length.
- 2. Activate both pieces of material by immersion in water bath 65°C (160°F).
- 3. Place the narrow Orficast More directly in the middle of the larger rectangle.
- 4. Briefly pat dry and bond the two layers together well.

Fabricate



- 1. Briefly pat dry and bond the two layers 4. With the patient standing in elbow together well.
- 2. Cut a small hole for the thumb in the center of one end, about 3-4 cm (1-1.5") from the border.
- 3. Opposite to this opening, cut an opening for the elbow's cubital fossa.



flexion with a forearm in neutral, place the material over the patient's thumb and pinch lightly near the fifth metacarpal.




5. Pull the material proximally and pinch together 6. Cut two small at the posterior elbow. Pinch lightly along the entire ulnar border of the forearm so that the material conforms to the shape. Mold carefully along the forearm structures.



pieces of material, activate in the splint pan and use these to reinforce the material directly over the lateral and medial epicondyles.

Position • The forearm is typically placed in neutral with the thumb pointing up - but the forearm can also be in full pronation or supination • Pearl: Stretch the material for maximum conformity. Pinch together lightly to allow for un pinching as soon as the shape has been captured.



7. Let the material harden and mark the trim lines. Undo all of the pinches and remove from the patient. Complete the trimming and apply straps.



Stay Connected

orfit EDGING STRIP

SOKT95K56

an de bordure

The community of Orfit users spans across the world and they love to share their most interesting patient cases and creative orthotic fabrication ideas.

Follow Orfit on Instagram for quick and interesting product updates. https://www.instagram.com/orfitindustries/

rehabilitation/

Learn and earn CEU's at the Orfit Academy

Learn from Orfit's Educational and Product Specialist, Debby Schwartz, OTD, OTR/L, CHT. Anytime, anywhere. Courses are free of charge. Payment is only required when taking the exam for accreditation.



www.orfit.com/academy

Where to buy



Follow Orfit on LinkedIn for quick and interesting product updates. https://www.linkedin.com/showcase/orfit-splinting-&-physical-

Join us now, learn from each other and support your colleagues. www.facebook.com/groups/Orfit.splinting

Register now to try one of our free courses to improve your orthotic fabrication skills.

> Scan the QR code or visit https://www.orfit.com/contact/physical-rehabilitation to learn where you can purchase Orfit products.

Visit WWW.ORFIT.COM for more information.

ORFIT INDUSTRIES

Vosveld 9A B-2110 Wijnegem | Belgium T (+32) (0)3 326 20 26 welcome@orfit.com

© in J f

ORFIT INDUSTRIES AMERICA

810 Ford Drive Norfolk, VA 23523 | USA T 516-935-8500 welcome@orfit.com



www.orfit.com