
3PK -3 Panel Knee Immobiliser

EN	KNEE ORTHOSIS
DE	KNIEORTHESE
FR	ORTHÈSE DU GENOU
ES	ORTESIS DE RODILLA
IT	ORTESI PER GINOCCHIO
NL	KNIEORTHESE
DK	KNÆORTOSE
FI	POLVIORTOOSI
SE	KNÄORTOS
NO	ORTOSE, KNE, ORDINÆR
PL	ORTEZA STAWU KOLANOWEGO
GR	ΝΑΡΘΗΚΑΣ ΓΟΝΑΤΟΣ

ENGLISH

Product Description

A three panel knee immobilisation orthosis with straps and stays to help immobilise the knee joint.

Usage Information (See diagram)

- 1) Remove the outer packaging and open all the straps on the knee immobiliser.
- 2) Position side panels to accommodate patient's leg size. (See Fig. 1)
- 3) Lay the knee immobiliser flat with the buckles and straps facing downwards. (See Fig. 2)
- 4) Immobilisation stays can be reshaped if required.
- 5) Slide the brace under the patient's leg. (See Fig. 3)
- 6) Fasten the straps to secure the brace to patient's leg. (See Fig. 4)

Indications for Use

- Conservative treatment for knee injuries.
- Maintain knee extension following surgery.
- Maintain knee extension following Trauma.
- Maintain knee extension following soft tissue injuries.

Contraindications

- Thrombosis of the deep veins of the leg and thigh.
- Do not apply in direct contact with broken skin.
- Any known allergy and/or hypersensitivity to any material listed in product composition.

Warnings and Precautions

We recommend the initial fitting of this device be conducted by a suitably qualified healthcare professional who will advise the period of use.

- Carefully read all instructions and warnings prior to use.
- Follow all instructions to ensure proper performance of the device.
- Do not use if liniments, ointments, gels, creams or any other substances have been applied to the affected area.
- Do not re-use for another patient, doing so risks cross infection and can compromise product integrity.
- Perform regular skin and circulation checks, especially patients with diabetes, vascular deficiencies and neurological conditions.
- Should any adverse reactions occur, please cease use and contact your healthcare professional or provider of this product.
- The durability of the device may be compromised by certain factors, e.g., objects with sharp edges or damage to the hook and loop fasteners. To prevent this the hook and loop fasteners should always be fastened when the device is not being worn or when being washed.

Washing and Care Instructions

1. Remove metal bars.
2. Close over all hook/loop straps.
3. Hand wash brace in cool water using mild detergent.
4. Rinse and dry flat.
5. Do not spin or tumble dry.
6. Do not use fabric conditioner at any time.
7. Before re-use, replace metal bars.
8. Check for any degradation that may affect use.
9. Periodically check hook/loop fasteners are operating securely and that the metal splint is secure in place.

Product Family Composition

Nylon, Cotton, Aluminium, Polyamide, Polyester, Polyurethane, POM, Steel.

Storage and Transport Conditions

Store in a cool, dry place out of direct sunlight, in the original packaging.

Recycling and Disposal

Packaging and constituent parts should be recycled or safely disposed in accordance with local or national laws.

Serious Incident

Report any serious incident to the manufacturer and the competent authority of the EU Member state, or country in which you reside.

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