
LSS - Lancaster Shoulder Sling

EN	SHOULDER ORTHOSIS
DE	SCHULTERORTHESE
FR	ORTHÈSE DE L'ÉPAULE
ES	ORTESIS DE HOMBRO
IT	ORTESI PER SPALLA
NL	SCHOUDERORTHESE
DK	SKULDERORTOSE
FI	OLKANIVELORTOOSI
SE	AXELORTOS
NO	ORTOSE, SKULDER
PL	ORTEZA STAWU BARKOWEGO
GR	ΝΑΡΘΗΚΑΣ ΩΜΟΥ

ENGLISH

Product Description

A universal shoulder orthosis to provide immobilisation holding the arm close to the body.

Usage Information (See diagram)

1. Place arm in sling and ensure the elbow is held securely in place.
2. Close the upper part of sling by attaching the closure strap (with no buckles on it) mid-way between wrist and elbow crease.
3. Attach the strap with the buckles on it, slightly above wrist as shown. Ensure buckle A is positioned at the back, facing the body. Ensure buckle B is at the top, facing upwards.
4. Take the shoulder strap up across the back.
5. Place the strap over the shoulder opposite the injured arm.
6. Fasten buckle B. Adjust the strap and secure with the hook/loop fastening. Excess strap length can be trimmed if required. (See step 10)
7. Take the body strap and attach it to the front of the arm sling as shown, using hook/loop fastening.
8. Bring the body strap around the back of the waist.
9. Fasten buckle A. Adjust the strap and secure with the hook/loop fastening. Excess strap length can be trimmed if required. (See step 10).
10. If shoulder and/or body strap length is too long, remove "Y" tabs and trim straps as necessary and reattach "Y" tab.
11. If required, adjust sling length by folding back using the lines marked on the sling. The folds can be kept in place using the circular hook pads provided.
12. Front and rear view of the completed shoulder brace fitting.

Indications for Use

- Can be used for an array of injuries and conditions where the restriction of shoulder and / or elbow movement is required such as:

- Dislocation.
- Strains and sprains.
- Pre and post operatively.
- Rehab stages.
- Other injuries.
- Brachial plexus or peripheral nerve.

Contraindications

- Do not apply the product in direct contact with broken skin.

Warnings and Precautions

We recommend the initial fitting of this device be conducted by a suitably qualified healthcare professional who will advise the period of use.

- Carefully read all instructions and warnings prior to use.
- Follow all instructions to ensure proper performance of the device.
- Do not use if liniments, ointments, gels, creams or any other substances have been applied to the affected area.
- Do not re-use for another patient, doing so risks cross infection and can compromise product integrity.
- Perform regular skin and circulation checks, especially patients with diabetes, vascular deficiencies and neurological conditions.
- Should any adverse reactions occur, please cease use and contact your healthcare professional or provider of this product.
- The durability of the device may be compromised by certain factors, e.g., objects with sharp edges or damage to the hook and loop fasteners. To prevent this the hook and loop fasteners should always be fastened when the device is not being worn or when being washed.

Washing and Care Instructions

1. Close over all hook/loop straps and fasten buckles and fasten buckles to the end.
2. Hand wash sling in cool water using mild detergent.
3. Rinse and dry flat.
4. Do not spin or tumble dry.
5. Do not use fabric conditioner at any time.
6. Check for any degradation that may affect use.
7. Periodically check hook/loop fasteners are operating securely.

Product Family Composition

Cotton, Nickel, Nylon, Polyester, Polypropylene, Polyurethane, POM, Steel.

Storage and Transport Conditions

Store in a cool, dry place out of direct sunlight, in the original packaging.

Recycling and Disposal

Packaging and constituent parts should be recycled or safely disposed in accordance with local or national laws.

Serious Incident

Report any serious incident to the manufacturer and the competent authority of the EU Member state, or country in which you reside.

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