
PAV - Pavlik Harness

EN	INFANT HIP HARNESS
DE	HÜFTGURT FÜR SÄUGLINGE
FR	HARNAIS DE HANCHE POUR ENFANT
ES	ARNÉS DE CADERA PARA BEBÉS
IT	ORTESI DELL'ANCA PER INFANZIA
NL	HEUPHARNAS VOOR BABY'S
DK	HOFTESELE TIL SPÆDBØRN
FI	VAUVAN LONKKAVALJAAT
SE	HÖFTSELE FÖR SPÄDBARN
NO	HOFTESELE FOR SPEDBARN
PL	ORTEZA BIODROWA DZIECIĘCA
GR	ΠΑΙΔΙΚΗ ΖΩΝΗ ΙΣΧΙΟΥ

ENGLISH

Product Description

To control hip flexion and abduction in the primary treatment of congenital dislocation of the hips.

Usage Information (See diagram)

- Open the colour-coded leg straps, shoulder straps unfasten the chest band and lay out the harness on a flat safe surface.
- Place the child’s back on the V portion of the harness.
- Fasten the chest band, just below the armpits (leaving a little clearance).
- Connect the shoulder straps to the front of the chest band.
- Place the child’s feet into the booties, ensuring the heel is located into the back of the bootie.
- Locate the leg straps through their buckles and pull through to desired position of flexion / abduction.
- Colour coding on the upper and lower straps prevents misapplication when re-fitting the harness.

Indications for Use

- Immobilisation following suspected congenital dislocation of the hip.
- Hip dysplasia.

Contraindications

- Do not use over open wounds.
- Any known allergy and/or hypersensitivity to any material listed in product composition.

Warnings and Precautions

We recommend the initial fitting of this device be conducted by a suitably qualified healthcare professional who will advise the period of use.

- Carefully read all instructions and warnings prior to use.
- Follow all instructions to ensure proper performance of the device.
- Do not use if liniments, ointments, gels, creams or any other substances have been applied to the affected area.
- Do not re-use for another patient, doing so risks cross infection and can compromise product integrity.
- Perform regular skin and circulation checks, especially patients with diabetes, vascular deficiencies and neurological conditions.
- Should any adverse reactions occur, please cease use and contact your healthcare professional or provider of this product.
- The durability of the device may be compromised by certain factors, e.g., objects with sharp edges or damage to the hook and loop fasteners. To prevent this the hook and loop fasteners should always be fastened when the device is not being worn or when being washed.

Washing and Care Instructions

- Close over all hook/loop straps.
- Hand wash in cool water using mild detergent.
- Rinse and dry flat.
- Do not spin or tumble dry.
- Do not use fabric conditioner at any time.
- Check for any degradation that may affect use.
- Periodically check hook/loop fasteners are operating securely.

Product Composition

Cotton, Elastane, Nylon, Polyamide, Polyester.

Storage and Transport Conditions

Store in a cool, dry place out of direct sunlight, in the original packaging.

Recycling and Disposal

Packaging and constituent parts should be recycled or safely disposed in accordance with local or national laws.

Serious Incident

Report any serious incident to the manufacturer and the competent authority of the EU Member state, or country in which you reside.

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