
PBS - Proflex Back Support

EN	LUMBAR SPINE ORTHOSIS
DE	LWS-ORTHESE
FR	ORTHÈSE DE LA COLONNE LOMBAIRE
ES	ORTESIS LUMBAR
IT	ORTESIS LUMBAR
NL	LUMBALE ORTHESE
DK	LÆNDEHVIRVELORTOSE
FI	LANNERANGAN ORTOOSI
SE	LÄNDRYGGSORTOS
NO	ORTOSE, LUMBAL
PL	ORTEZA KRĘGOSŁUPA ŁĘDŹWIOWEGO
GR	ΝΑΡΘΗΚΑΣ ΟΣΦΥΙΚΗΣ ΜΟΙΡΑΣ ΤΗΣ ΣΠΟΝΔΥΛΙΚΗΣ ΣΤΗΛΗΣ

ENGLISH

Product Description

Soft-good/semi-rigid designed to restrict motion through elastic or semi-rigid construction. May be suitable for the symptomatic relief of acute and chronic low back pain.

Usage Information (See diagram)

- Wrap the binder around your front at the level of your lower rib cage, just below your chest.
- Using your left hand, stretch the end of the binder across your body and hold against your ribs.
- Using your right hand, stretch the fastener end of the binder and across your body.
- Secure the fastener strip on padding at your desired level of compression and support.

Indications for Use

• Acute or chronic lower back pain.

Contraindications

• Do not use over open wounds.
• Any known allergy and/or hypersensitivity to any material listed in product composition.

Warnings and Precautions

We recommend the initial fitting of this device be conducted by a suitably qualified healthcare professional who will advise the period of use.
• Carefully read all instructions and warnings prior to use.
• Follow all instructions to ensure proper performance of the device.
• Do not use if liniments, ointments, gels, creams or any other substances have been applied to the affected area.
• Do not re-use for another patient, doing so risks cross infection and can compromise product integrity.
• Perform regular skin and circulation checks, especially patients with diabetes, vascular deficiencies and neurological conditions.
• Should any adverse reactions occur, please cease use and contact your healthcare professional or provider of this product.
• The durability of the device may be compromised by certain factors, e.g., objects with sharp edges or damage to the hook and loop fasteners. To prevent this the hook and loop fasteners should always be fastened when the device is not being worn or when being washed.

Washing and Care Instructions

- Close over all hook/loop straps.
- Hand wash in cool water using mild detergent.
- Rinse and dry flat.
- Do not spin or tumble dry.
- Do not use fabric conditioner at any time.
- Check for any degradation that may affect use.
- Periodically check hook/loop fasteners are operating securely.

Product Family Composition

Cotton, Elastane, Natural Rubber, Nylon, Polyamide, Polyester, Polyethylene, Steel, Thermoplastic, Viscose.

Storage and Transport Conditions

Store in a cool, dry place out of direct sunlight, in the original packaging.

Recycling and Disposal

Packaging and constituent parts should be recycled or safely disposed in accordance with local or national laws.

Serious Incident

Report any serious incident to the manufacturer and the competent authority of the EU Member state, or country in which you reside.

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©2023 Promedics Orthopaedics Ltd. All Rights Reserved.