





SI - Shoulder Immobiliser

EN	SHOULDER ORTHOSIS
DE	SCHULTERORTHESE
FR	ORTHÈSE DE L'ÉPAULE
ES	ORTESIS DE HOMBRO
IT	ORTESI PER SPALLA
NL	SCHOUDERORTHESE
DK	SKULDERORTOSE
FI	OLKANIVELORTOOSI
SE	AXELORTOS
NO	ORTOSE, SKULDER
PL	ORTEZA STAWU BARKOWEGO
GR	ΝΑΡΘΗΚΑΣ ΟΜΟΥ

ENGLISH

Product Description

A universal low-profile shoulder orthosis to provide shoulder immobilisation.

Usage Information (See diagram)

- Position body wrap around waist and secure with hook/loop fastening.
- Position upper arm wrap above elbow joint and secure with hook/loop fastening. Once secure use hook/loop fastening to attach to body wrap.
- Secure wrist wrap around wrist with hook/loop fastening ensuring buckle is facing upwards.
- Secure shoulder strap to rear of body wrap before gently pulling strap up over unaffected shoulder and locking the buckle (use scissors to trim shoulder strap for comfortable fit).
- Position thumb wrap around thumb, securing with hook/loop fastening below buckle.
- Take care not to over tighten straps, re-adjust for comfort as required while ensuring orthosis is in correct position.

Indications for Use

Can be used for an array of injuries and conditions where the restriction of shoulder movement is required such as:

- Dislocation.
- Strains and sprains.
- Pre and post operatively.
- Rehab stages.
- Post fracture management.
- Other injuries.

It is recommended that individuals seek guidance from a medial practitioner regarding the most appropriate orthosis.

Contraindications

- Do not use over open wounds.
- Any known allergy and/or hypersensitivity to any material listed in product composition.
- Caution should be used if prescribed to patients with diabetes, vascular deficiency and neuropathy.

Warnings and Precautions

We recommend the initial fitting of this device be conducted by a suitably qualified healthcare professional who will advise the period of use.

- Carefully read all instructions and warnings prior to use.
- Follow all instructions to ensure proper performance of the device.
- Do not use if liniments, ointments, gels, creams or any other substances have been applied to the affected area.
- Do not re-use for another patient, doing so risks cross infection and can compromise product integrity.
- Perform regular skin and circulation checks, especially patients with diabetes, vascular deficiencies and neurological conditions.
- Should any adverse reactions occur, please cease use and contact your healthcare professional or provider of this product.
- The durability of the device may be compromised by certain factors, e.g., objects with sharp edges or damage to the hook and loop fasteners. To prevent this the hook and loop fasteners should always be fastened when the device is not being worn or when being washed.

Washing and Care Instructions

- Close over all hook/loop straps.
- Hand wash in cool water using mild detergent.
- Rinse and dry flat.
- Do not spin or tumble dry.
- Do not use fabric conditioner at any time.
- Check for any degradation that may affect use.
- Periodically check hook/loop fasteners are operating securely.

Product Family Composition

Elastane, Nylon, Polyamide, Polyester, Polyurethane, POM, Steel, Spandex.

Storage and Transport Conditions

Store in a cool, dry place out of direct sunlight, in the original packaging.

Recycling and Disposal

Packaging and constituent parts should be recycled or safely disposed in accordance with local or national laws.

Serious Incident

Report any serious incident to the manufacturer and the competent authority of the EU Member state, or country in which you reside.

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