UFO - Plantar Fasciitis Orthosis

EN ANKLE/FOOT ORTHOSIS

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ES ORTESIS DE TOBILLO/PIE

IT ORTESI DI CAVIGLIA/PIEDE

NL ENKEL-/VOETORTHESE

DK ANKEL-/FODORTOSE

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
SE ORTOS, FOT/FOTLED

NO ORTOSE, ANKEL/FOT

PL ORTEZA GOLEŃ/STOPA

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Ankle orthosis night splint for plantar fasciitis.

Usage Information (See diagram)

- Open straps 1 through to 3 and place foot into the boot.
- Fasten strap 1.
- Tighten strap 1.
- Fasten and tighten all remaining straps.
- Adjust dorsiflexion strap as required. (Figure 6)

Indications for Use

- Alleviate pain associated with plantar fasciitis.
- Treatment of Achilles tendonitis.

Contraindications

- Do not attempt to walk whilst using this device.
- Not to used directly over open wounds.

Warnings and Precautions

We recommend the initial fitting of this device be conducted by a suitably qualified healthcare professional who will advise the period of use.

- Carefully read all instructions and warnings prior to use.
- Follow all instructions to ensure proper performance of the device.
- Do not use if liniments, ointments, gels, creams or any other substances have been applied to the affected area.
- Do not re-use for another patient, doing so risks cross infection and can compromise product integrity.
- Perform regular skin and circulation checks, especially patients with diabetes, vascular deficiencies and neurological conditions.
- Should any adverse reactions occur, please cease use and contact your healthcare professional or provider of this product.
- The durability of the device may be compromised by certain factors, e.g., objects with sharp edges or damage to the hook and loop fasteners. To prevent this the hook and loop fasteners should always be fastened when the device is not being worn or when being washed.

Washing and Care Instructions

- Close over all straps and remove liner.
- Hand wash liner in cool water using mild detergent.
- Rinse and dry flat.
- Do not spin or tumble dry.
- Do not use fabric conditioner at any time.
- Wipe down boot with cloth or cleaning wipe.
- Check for any degradation that may affect use.
- Periodically check straps are operating securely.

Product Family Composition

EVA, Polyester, Polypropylene, Polyurethane, POM.

Storage and Transport Conditions

Store in a cool, dry place out of direct sunlight, in the original packaging.

Recycling and Disposal

Packaging and constituent parts should be recycled or safely disposed in accordance with local or national laws.

Serious Incident

Report any serious incident to the manufacturer and the competent authority of the EU Member state, or country in which you reside.

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©2023 Promedics Orthopaedics Ltd. All Rights Reserved.