

Product Details

This premium knee wrap is ideal for patients with mild-to-moderate Osteoarthritis. The wraparound design features a breathable, high compression material with patient friendly closure tabs for easy application.

+ **Breathable**

High compression, breathable wrap with patient friendly closure tabs for easy and correct application.

+ **Range of Motion**

ROM settings can be adjusted using inserts supplied.

+ **Adjustable**

Double adjustable hinge allows for valgus or varus correction. Various settings above and below the knee joint allows improved fit, allows for a more even distribution.

+ **Effective Off-Loading**

Two point hinge adjustment above and below the knee creates effective valgus or varus off-loading.



Scan QR code to
watch fitting video
www.promedics.co.uk

Tel: 01475 746 400

Email: enquiries@promedics.co.uk

www.promedics.co.uk

Promedics Orthopaedics Ltd
Port Glasgow Industrial Estate
Gareloch Road, Port Glasgow
Scotland, PA14 5XH



Life Changing Innovation

 **promedics**
an Aspen company

OA knee+

NEW



 **promedics**
an Aspen company

OA Knee+



Size	Code
Left Medial	PROAP/L*
Right Medial	PROAP/R*

*For lateral OA, use /L on the right leg and /R on the left leg.



OA Wrap



Product Details

This premium brace is ideal for patients with mild-to-moderate Osteoarthritis. The low profile, universal design allows for less inventory and better patient compliance.

- + Comfort**
Comfortable, non-slip pads to help prevent migration.
- + Range of Motion**
Adjustable stops allow you to set the range of motion for appropriate treatment during continuum of care
- + Two Point Hinge Adjustment**
Hinge adjustment for even pressure distribution and comfort.
- + One-Step Closure**
Proprietary one-step closure system for easy donning and doffing.

Size	Code Left Medial*	Code Right Medial*	Thigh (cm)	Knee (cm)	Calf (cm)
S	PROAP/13	PROAP/23	39 - 47	33 - 35	30.5 - 35.5
M	PROAP/15	PROAP/25	47 - 53.5	35.5 - 38	35.5 - 40.5
L	PROAP/17	PROAP/27	53.5 - 59.5	38 - 43	40.5 - 47
XL	PROAP/18	PROAP/28	59.5 - 67.5	43 - 48.5	47 - 50.5
XXL	PROAP/19	PROAP/29	48.5 - 53.5	48.5 - 53.5	50.75 - 56

*For lateral OA, use corresponding left code on the right leg and corresponding right code on the left leg.