

Carbon AFO Medial

Product Code: KTAFO1/M*

The Medial Model is suited for use on patients with a supinated ankle (leaning outward) as the strut is away from the ankle moving laterally.

- + Lightweight & Strong
- + Strapping Fitted & Fully Trimmable
- + Removable / Replaceable Ankle Pads
- + Ideal Stem Placement
- + Assists Foot Drop

Size	Left	Right
S	KTAFO1/M/SL	KTAFO1/M/SR
M	KTAFO1/M/ML	KTAFO1/M/MR
L	KTAFO1/M/LL	KTAFO1/M/LR
XL	KTAFO1/M/XL	KTAFO1/M/XLR



Carbon AFO

Lateral

Product Code: KTAFO1/L*

The Lateral Model is suited for use on patients with a pronated ankle (leaning inward) as the strut is away from the ankle moving medially.

- + Lightweight & Strong
- + Strapping Fitted & Fully Trimmable
- + Removable / Replaceable Ankle Pads
- + Ideal Stem Placement
- + Assists Foot Drop

Size	Left	Right
S	KTAFO1/L/SL	KTAFO1/L/SR
M	KTAFO1/L/ML	KTAFO1/L/MR
L	KTAFO1/L/LL	KTAFO1/L/LR
XL	KTAFO1/L/XL	KTAFO1/L/XLR

