

Oedema Gloves

Made from a soft, silky nylon lycra with seams only on the outside, our Oedema Gloves are available in two wrist lengths and two finger styles. Provides warmth and compression for soothing soft tissue oedema and arthritic joint pain. Features an open slit in the cuff for easy application and non-constricting fit.

Provides 15 - 25mmHg of gentle and even compression to oedematous tissue. Wrist length measures 25mm proximal to wrist crease. Over the wrist length measures 76mm proximal to wrist crease.



Full Finger Code (Left)	Full Finger Code (Right)	Full Finger Over Wrist Code (Left)	Full Finger Over Wrist Code (Right)	Open Finger Over Wrist Code (Left)	Open Finger Over Wrist Code (Right)	Size
JOG/1/XXSL	JOG/1/XXSR	JOG/2/XXSL	JOG/2/XXSR	JOG/3/XXSL	JOG/3/XXSR	XXS
JOG/1/XSL	JOG/1/XSR	JOG/2/XSL	JOG/2/XSR	JOG/3/XSL	JOG/3/XSR	XS
JOG/1/SL	JOG/1/SR	JOG/2/SL	JOG/2/SR	JOG/3/SL	JOG/3/SR	S
JOG/1/ML	JOG/1/MR	JOG/2/ML	JOG/2/MR	JOG/3/ML	JOG/3/MR	M
JOG/1/LR	JOG/1/LR	JOG/2/LL	JOG/2/LR	JOG/3/LL	JOG/3/LR	L