

# Neoprene Thumb Spica

A Neoprene thumb spica featuring a removable flexible stay which restricts movement at the CMC and MCP joints, leaving the IP joint free.



Combining the advantages of a thumb splint with a wrist splint, the Thumb Spica provides great finger support to relieve pain associated with a number of different issues.

This is a universal support, meaning it can be worn on either hand. This is particularly useful if you suffer from a recurring condition, allowing you to wear the thumb spica on either hand requiring support.

Code (Left)	Code (Right)	Size
NTS/13	NTS/23	S
NTS/15	NTS/25	M
NTS/17	NTS/27	L
NTS/18	NTS/28	XL