

Ventilated Thumb Restriction Splint

A ventilated, low profile wrap around support providing adjustable support for the thumb CMC joint in the treatment of osteoarthritis and medial subluxation. The shaped strap wraps around the thumb CMC joint and the first metacarpal for targeted CMC joint support.

An extensive size range covers paediatric to adult sizes and plus sizes are a full size bigger in thumb circumference to fit oedematous thumb conditions. Function and fit without bulk.

Provides wrist and thumb support for an array of injuries, conditions, pre-and post-operations such as:

- + CMC joint arthritis
- + CMC joint instability
- + Tendinitis





Designed with a ventilated, comfortable material and a simple hook-and-loop closure. The low-profile wrap around the thumb CMC joint and first metacarpal provides direct support of the CMC joint.

Code (Left)	Code (Right)	Size
TRSV/11	TRSV/21	XXS
TRSV/12	TRSV/22	XS
TRSV/13	TRSV/23	S
TRSV/13+	TRSV/23+	S+
TRSV/15	TRSV/25	M
TRSV/15+	TRSV/25+	M+
TRSV/17	TRSV/27	L
TRSV/17+	TRSV/27+	L+
TRSV/18	TRSV/28	XL