

Orfilight Thumb Post Splint

This orthosis gives an extremely comfortable immobilisation and stabilisation of the complete thumb column. It can serve for the treatment of De Quervain tenosynovitis, rheumatoid arthritis, basal joint osteoarthritis, thumb sprains and pre- or post-operative thumb management.

Weighing 15-20% less than other splinting materials, Orfilight[®] has excellent elasticity, significant (yet slow) memory and a soft feel. Offering high stretch, it allows for intimate moulding to anatomical contours. Self-adhesive for reliable bonding and easy fixation of accessories. Orfilight[®] has maximum stretch without tearing and is completely remouldable.



Size	Beige	Black	Material & Perforation
S	-	35810Z/L/NS	1.6mm Micro Perfo
M	35811/L	35811Z/L/NS	1.6mm Micro Perfo
L	35812/L	35812Z/L/NS	1.6mm Micro Perfo