

Horizon Pro Lower Spine Braces

Backed by clinical research and 10 years of trusted technology, the Horizon PRO provides an updated approach to pain management. Its inelastic design combined with targeted compression offers immediate relief from symptoms caused by injury, muscle fatigue, and spasms.

One Size Adjustable

Comfortably fits waists ranging from 61-127cm and up to 178cm with the extension panels. Simple pull-through mechanism offers multiple sizes in one brace, enabling each brace to comfortably fit various patient shapes and sizes.



Horizon PRO 464

Horizon PRO 456



Manage symptoms of low back pain and add trunk stability with market-leading technology

Integrated Anterior Panels

Anterior panels overlap to provide enhanced midline support. Featuring Aspen's innovative FlexTabs™, these panels disperse pressure and naturally conform to the anatomy for added comfort.

Low Profile Design with Reliable Structural Integrity

The Horizon PRO 639 features a tall, rigid anterior panel and overlapping lateral/anterior panels that provide flexion and lateral motion restriction without compromising comfort. The xiphoid extension can be angled as needed or removed for shorter patients.

Ergonomic Pull Tabs

Comfortable and easy to use for patients to achieve their desired level of compression.

Easy Access to Panels

Improve customisation in the clinic.



Custom Support

When tightening, the firm but flexible back panel naturally conforms to the patient's lordosis to provide added structural support.

Integrated Pockets

Provides easier application and ensures a snug fit.

Height Adjusts for Every Body

Telescoping shoulder straps on the Horizon PRO 456 provide 33cm of height adjustment to fit varying patient anatomies. In addition, circumferential adjustability enables a customised fit for every body.

Independent Tightening Mechanism

Aspen's patented 4:1 mechanically advantaged SlickTrack™ System features an independent upper and lower tightening mechanism, providing direct compression to support tight muscles.



The Right Brace for Multiple Needs

Regardless of which configuration a patient requires initially, all Horizon PRO braces can be stepped down to lowest-profile option for the treatment of reoccurring ailments.

FITS WAIST CIRCUMFERENCE RANGING FROM 61-127CM AND UP TO 78CM WITH HORIZON EXTENSION PANEL

Product	Code
Horizon™ PRO 464 TLSO	994760
Horizon™ PRO 456 TLSO	994740
Horizon™ PRO 639 LSO	994735
Horizon™ PRO 637 LSO	994730

Product	Code
Horizon™ PRO 631 LSO	994720
Horizon™ PRO 627 Lumbar	994710
Horizon™ Extension Panels	993745
Therapy Pack	990005