

AirX Thumb Restriction Splint

A low profile AirX wrap that allows uninterrupted finger function and provides adjustable support to the thumb CMC joint. Excellent in the treatment of osteoarthritis and medial subluxation. The adjustable strap can be wrapped around the first CMC and MCP joints to provide targeted joint support.

Provides wrist and thumb support for an array of conditions such as:

- + CMC joint arthritis
- + CMC joint instability
- + Tendonitis



Size	Left	Right	Wrist Circumference
Infant	ATRS/11	ATRS/21	11cm - 14cm
Child	ATRS/12	ATRS/22	14cm - 17cm
S	ATRS/13	ATRS/23	15cm - 18cm
S+	ATRS/13+	ATRS/23+	17cm - 18cm
M	ATRS/15	ATRS/25	18cm - 20cm
M+	ATRS/15+	ATRS/25+	20cm - 21cm
L	ATRS/17	ATRS/17	20cm - 23cm
L+	ATRS/17+	ATRS/17+	21cm - 24cm
XL	ATRS/18	ATRS/28	23cm - 25cm