

Abducted Thumb Wrap

The Abducted Thumb Wrap is an ideal choice as the closed thumb support holds the joint in the natural abduction position to relieve strain. The removable stay helps immobilise the palm and thumb, while the neoprene material ensures maximum comfort while wearing the brace.

Provides wrist and thumb support for an array of injuries, conditions, pre-and post-operations such as:

- + Arthritis
- + Sprains and strains
- + Soft tissue injuries and conditions
- + Instability



Easy hook / loop closure for simple application and removal



Palmar and thumb extension bar shaped for best fit, providing additional support for the first CMC and MCP joints

Size	Left	Right	Wrist Circumference
XS	ATW/12	ATW/22	11.5cm - 13cm
S	ATW/13	ATW/23	14cm - 15cm
M	ATW/15	ATW/25	16.5cm - 18cm
L	ATW/17	ATW/27	20cm - 21.5cm
XL	ATW/18	ATW/28	23cm - 24cm