

Single Armchair Finger Splint



A dynamic splint with coil design for PIP and DIP joint flexion contractures and Boutonniere's deformity. Maintains assisted extension at the PIP joint while allowing active flexion and can be used for soft tissue injury affecting the PIP joint with a possible flexion deformity.

Provides finger support / extension for an array of injuries, conditions, pre and post operations such as:

- + PIP & DIP joint flexion contractures
- + Extensor tendon injuries / conditions (Boutonniere's deformity)
- + Soft tissue conditions & injuries



Size	Code	Circumference at PIP
XS	CAP/2XS	4cm - 5cm
S	CAP/2S	5.5cm - 6cm
M	CAP/2M	6cm - 6.5cm
L	CAP/2L	6.5cm - 7cm
XL	CAP/2XL	7cm - 7.5cm