

Promedics Cryo Cold Therapy System

5X
WRAP OPTIONS
AVAILABLE

The Promedics Cryo Cold Therapy System is an easy to use device that uses gravity to circulate ice cold water providing cold therapy in the treatment of injuries. Using thermal insulation enhanced technology, it reduces cooling loss and prevents heat entering the system, staying colder for longer. **All five wraps are universal, one size fits all.**

Promedics Cryo is supplied with tube included. It is ideal for post-operative recovery, trauma, swelling and rehabilitation.



SHOULDER WRAP



ELBOW WRAP



KNEE WRAP



ANKLE WRAP



BACK WRAP



TUBE

Product	Code	NHSSC Code	Size	Weight
Cooler	CRYO/COOLER	FBG85004	320mm x 185mm (diameter)	1077g
Shoulder Wrap	CRYO/SW/10	FBG85001	965mm x 465mm	255g
Elbow Wrap	CRYO/EW/10	FBG85005	610mm x 385mm	191g
Knee Wrap	CRYO/KW/10	FBG85002	835mm x 345mm	234g
Ankle Wrap	CRYO/AW/10	FBG85000	725mm x 375mm	267g
Back Wrap	CRYO/BW/10	FBG85003	1200mm x 305mm	284g
Tube	CRYO/TUBE	FBG85066	110mm x 17mm (diameter)	80g

CRYOPRO/SSB - 5/24