

Deltaform Plus Wrist Brace

The Deltaform Plus Wrist Brace is a 23cm long woven elastic/cotton support, with an aluminium stay inside a vinyl pocket and a longitudinal seam for excellent conformity to the wrist and forearm. Extends to two thirds of the forearm for radial/ulnar stability. **Available in beige and black.**

It offers support, compression and stabilisation. The gentle compression provided by the elastic helps to reduce swelling and assists with the natural healing process. The simple hook loop fastening, which fastens securely to provide a personalised fit and easy application.

Provides wrist and thumb support for an array of injuries, conditions, pre-and post-operations such as:

- + Arthritis
- + Sprains and strains
- + Soft tissue injuries and conditions
- + Instability
- + Post fracture management



Size	Black - Left	Black - Right	Beige - Left	Beige - Right	Wrist Circ.
XS	DFP/12/BK	DFP/22/BK	DFP/12	DFP/22	11.5cm - 13cm
S	DFP/13/BK	DFP/23/BK	DFP/13	DFP/23	14cm - 15cm
M	DFP/15/BK	DFP/25/BK	DFP/15	DFP/25	16.5cm - 18cm
L	DFP/17/BK	DFP/27/BK	DFP/17	DFP/27	20cm - 21.5cm
XL	DFP/18/BK	DFP/28/BK	DFP/18	DFP/28	23cm - 24cm