

Evans

Abduction Sling

This is a low cost universal option that utilizes the popular Lancaster Shoulder Immobiliser. In addition to the comfortable sling with shoulder and body straps – the Evans Abduction Sling incorporates a 15° abduction pad that is easily secured to the inside of the sling. The inner profile of the abduction pad is shaped to rest comfortably against the patient.

- + Rotator cuff repairs
- + Glenohumeral dislocations / subluxation
- + Posterior / Anterior capsule repairs
- + Soft tissue repairs / strains



Fast release buckles for easy patient removal & re-application



Removable Hook/Loop 'Y' - Tabs



Adjust sling length using fold lines and circular hook pads

Size	Code
Universal	EAS