

# Elastic Thumb Spica

Providing protection, support and comfort while your thumb heals, the Elastic Thumb Spica is made from highly breathable fabric that prevents moisture build-up, leaving your hand feeling fresh all day. The elastic material imparts a gentle compression to the injured area, improving circulation and helping to reduce swelling.

It provides adjustable support to the thumb carpometacarpal (CMC) and metacarpal phalangeal (MCP) joint. A flexible removable stay restricts movement at the CMC and MCP joints, leaving the interphalangeal (IP) joint free.

Provides thumb CMC and MCP joint support for an array of injuries, conditions, pre- and post-operations such as:

- + Arthritis
- + Sprains and strains
- + Soft tissue injuries and conditions
- + Instability



Convenient hook / loop closures make it easy to adjust the support, and are ideal for users with limited dexterity.



Removable, flexible thumb stay restricts movement at the CMC and MCP joints for accelerated healing.

| Size | Left            | Right           | Wrist Circumference |
|------|-----------------|-----------------|---------------------|
| XS   | <b>ETSLF/12</b> | <b>ETSLF/22</b> | 11cm - 13cm         |
| S    | <b>ETSLF/13</b> | <b>ETSLF/23</b> | 14cm - 15cm         |
| M    | <b>ETSLF/15</b> | <b>ETSLF/25</b> | 15cm - 16.5cm       |
| L    | <b>ETSLF/17</b> | <b>ETSLF/27</b> | 16.5cm - 18cm       |
| XL   | <b>ETSLF/18</b> | <b>ETSLF/28</b> | 18cm - 19cm         |