

Elastic Thumb Spica - Modified

One-way stretch elastic supports the wrist and thumb, leaving the palm completely free. Angled hook/loop attachment at wrist and thumb ensures perfect fit. The 10cm removable, flexible thumb stay restricts movement at the CMC and MCP joints, leaving the IP joint free for an array of injuries and conditions such as;

- + Arthritis
- + Sprains and strains
- + Soft tissue injuries & conditions
- + Instability



Size	Left	Right	Wrist Circumference
XS	ETSMLF/12	ETSMLF/22	11cm - 13cm
S	ETSMLF/13	ETSMLF/23	14cm - 15cm
M	ETSMLF/15	ETSMLF/25	15cm - 16.5cm
L	ETSMLF/17	ETSMLF/27	16.5cm - 18cm
XL	ETSMLF/18	ETSMLF/28	18cm - 19cm