

FLX

Lower Spine Braces

Enhancing the proven design and function of the OTS product line, the FLX line of inelastic braces offers a conservative treatment solution for patients experiencing lower back pain. Structural improvements have been added to enhance stability and support while maintaining a sleek, breathable profile. These one-size adjustable braces are flexible by design, providing targeted compression and stability to support patients and aid in pain relief.



FLX 456 TLSO

FLX 464 TLSO

FLX 631 LSO



Integrated Anterior Panels

Anterior panels overlap to provide enhanced midline support. Featuring Aspen's innovative FlexTabs™, these panels disperse pressure and naturally conform to the anatomy for added comfort.

Improved Structural Integrity

Added reinforced materials to belt to improve support, stability, and compression, all while providing uncompromised comfort.


Added Access for Modification

Anterior and posterior panels are now fully accessible and removable to better accommodate a broader range of patient anatomies and allow for the utilisation of custom codes.



Custom Support

When tightening, the firm but flexible back panel naturally conforms to the patient's lordosis to provide added structural support.



Independent Tightening Mechanism

Aspen's patented 4:1 tightening mechanism provides targeted compression where needed to support tight muscles.



Adjusts for Every Body

Movable shoulder pads allow you to adjust the brace to your patient's specific needs. Optimising comfort and fit.



One Size Adjustable

Simple pull through mechanism offers multiple sizes in one brace, enabling each brace to comfortably fit various anatomies. This means the correct size is always on hand, minimizing inventory requirements and significantly reducing costs.



The Right Brace for Multiple Needs

Improvements to better accommodate your practice needs. Added posterior flap and anterior openings to allow easy access and removal of the panels, giving you better flexibility with fitment.

← FITS WAIST CIRCUMFERENCE RANGING FROM 71 - 137CM AND UP TO 162CM WITH FLX UNIVERSAL EXTENSIONS PANEL →

Product	Code
FLX 464 TLSO	514643
FLX 456 TLSO	514563
FLX 631 LSO	516313
FLX Universal Extension Panel	506505