

Jura Basic Mallet Finger Splint

A finger orthosis for immobilisation of the distal interphalangeal (DIP) joint of the digits in neutral. For digit injuries / conditions that require immobilisation / support of the distal interphalangeal (DIP) joint in neutral such as;

- + Mallet injuries
- + Fractures
- + Joint instabilities
- + Nail bed injuries



Size	Code	PIP Circumference
Small	JBPM/1	4.5cm - 5cm
Small/Medium	JBPM/2	5cm - 5.5cm
Medium	JBPM/3	6cm - 6.5cm
Medium/Large	JBPM/4	6.5cm - 7cm
Large	JBPM/5	7cm - 7.5cm
Large+	JBPM/5.5	7.5cm - 8cm
Large/X-Large	JBPM/6	8cm - 8.5cm
X-Large	JBPM/7	8.5cm - 9cm