

# Jura Clavicle Splint

A traditional figure-of-eight bandage design to provide support for immobilisation of the shoulder. The brace has soft padded shoulder straps with return loops and a padded back pad. Comfort pads can be positioned anywhere along the length of the straps; for example in the axilla area. Straps can be easily trimmed to length and refastened using the hook Y-tab closures.

- + For the conservative treatment of clavicle (collarbone) fractures / injury
- + For postoperative and post traumatic treatment
- + Shoulder retraction
- + Posture support



Size	Code
Universal	JCS/10