

Jura Universal Wrist Brace

Dual compression elastic for improved support and fit. A two piece brace with longitudinal seam conforms perfectly to the wrist, hand and forearm. Hard wearing webbing bar casing with aluminium palmar bar. Hook/loop fastening.

Provides wrist immobilisation / support for an array of injuries, conditions, pre- and post-operations such as:

- + Sprains and strains
- + Soft tissue injuries and conditions
- + Arthritis
- + Instability
- + Peripheral nerve injuries and conditions
- + Post fracture management



Size	Left	Right
Paediatric	JUWB/LP	JUWB/RP
Universal	JUWB/L	JUWB/R