

# Jura

## Wrist Thumb

Dual compression elastic for improved support and fit. A two-piece wrist thumb brace with longitudinal seam conforms perfectly to the wrist, hand and forearm. Aluminium palmar bar and thumb stay provides excellent support to wrist and first CMC & MCP joints.

The thumb stay can be easily adjusted for the required degree of abduction and the thumb closure is easily adjusted for a comfortable fit.

**Available in standard (18cm) or long (23cm) versions.**



Size	Left (Standard)	Right (Standard)	Left (Long)	Right (Long)	Wrist Circ.
XS	<b>JWT/12</b>	<b>JWT/22</b>	<b>JWTL/12</b>	<b>JWTL/22</b>	13cm - 14cm
S	<b>JWT/13</b>	<b>JWT/23</b>	<b>JWTL/13</b>	<b>JWTL/23</b>	14cm - 16cm
M	<b>JWT/15</b>	<b>JWT/25</b>	<b>JWTL/15</b>	<b>JWTL/25</b>	16cm - 19cm
L	<b>JWT/17</b>	<b>JWT/27</b>	<b>JWTL/17</b>	<b>JWTL/27</b>	19cm - 20cm
XL	<b>JWT/18</b>	<b>JWT/28</b>	<b>JWTL/18</b>	<b>JWTL/28</b>	20cm - 22cm