

Jura Wrist / Thumb Long

A wrist thumb orthosis with malleable aluminium thumb stay to hold the thumb firmly in abduction, the Jura Wrist / Thumb Brace supports the thumb MP joint and restricts movement at the CMC joint.

Provides wrist support for an array of injuries, conditions, pre-and post-operations such as:

- + Arthritis
- + Sprains and strains
- + Soft tissue injuries and conditions
- + Instability
- + Post fracture management



Size	Left	Right	Wrist Circ.
XS	JWTL/12	JWTL/22	< 15cm
S	JWTL/13	JWTL/23	15cm - 17cm
M	JWTL/15	JWTL/25	17cm - 19cm
L	JWTL/17	JWTL/27	19cm - 22cm
XL	JWTL/18	JWTL/28	23cm - 25cm