

# Jura Wrist / Thumb Brace

A wrist thumb orthosis with malleable aluminium thumb stay to hold the thumb firmly in abduction, supporting the thumb MP joint and restricting movement at the CMC joint.

Provides wrist support for an array of injuries, conditions, pre-and post-operations such as:

- + Arthritis
- + Sprains and strains
- + Soft tissue injuries and conditions
- + Instability
- + Post fracture management



| Size | Left    | Right   | Wrist Circumference |
|------|---------|---------|---------------------|
| S    | JWTO/13 | JWTO/23 | 12cm                |
| M    | JWTO/15 | JWTO/25 | 14cm                |
| L    | JWTO/17 | JWTO/27 | 16cm                |