

Lacer Wrist

Long / Standard

A premium orthosis for support and stabilisation of the wrist, with a conformable stay and wrist strap to aid stabilisation. The lacing closure allows for easy one-handed application, while the palmar pad provides additional comfort and support.

Provides wrist support for an array of injuries and conditions, such as;

- + Arthritis
- + Sprains, strains & instability
- + Soft tissue injuries & post fracture management



The lacer closure system allows for easy one-hand fastening with even tension.



The removable aluminium palmar bar is malleable for a custom, comfortable fit.



Available in both long and standard versions, with a range of sizes to suit all conditions.



The cushioned palmar pad provides additional comfort and support.

Size	Long - Left	Long - Right	Short - Left	Short - Right	Wrist Circumference
XXS	-	-	LCRS/11	LCRS/21	< 15cm
XS	-	-	LCRS/12	LCRS/22	15cm - 17cm
S	LCRL/13	LCRL/23	LCRS/13	LCRS/23	17cm - 20cm
M	LCRL/15	LCRL/25	LCRS/15	LCRS/25	20cm - 23cm
L	LCRL/17	LCRS/27	LCRS/17	LCRS/27	22cm - 27cm
XL	LCRL/18	LCRL/28	LCRS/18	LCRS/28	25cm - 31cm