

# Traverse Air Walker Tall & Short

The Traverse Air Walker is a circumferential walking boot that utilises dual air pumps for targeted, comfortable compression. Walker boots are designed to help prevent further damage and promote quicker healing of the injured site. The walker boot helps the user to become mobile again by protecting the site and slowly allowing weight to be put onto the injured foot/ankle.

- + Comfortable and lightweight
- + Easy to apply and remove
- + Air pump pressure can be controlled from either side of the ankle
- + Rocker base promotes a natural and stable gait
- + Shock absorbing heel design



Size	Men's Shoe Size	Women's Shoe Size
XS	Up to 3.5	Up to 4.5
S	3.5 - 6.5	2.5 - 5.5
M	6.5 - 9	5.5 - 8
L	9 - 11	8 - 10
XL	11 - 14.5	10 - 13.5

Traverse Air Walker - Short				
XS	S	M	L	XL
653001	653002	653003	653004	653005

Traverse Air Walker - Tall				
XS	S	M	L	XL
653011	653012	653013	653014	653015

Achilles Wedge		Extended Straps
Left	Right	Universal
JWW/RL	JWW/RR	JWA/EXS

Traverse Replacement Liner - Short				
XS	S	M	L	XL
653301	653302	653303	653304	653305

Traverse Replacement Liner - Tall				
XS	S	M	L	XL
653306	653307	653308	653309	653310

Traverse Hygiene Cover				
XS	S	M	L	XL
653401	653402	653403	653404	653405