

# LNK-Up

## NEW PRODUCT

The LNK-Up is designed for flaccid paralysis drop foot by raising the forefoot when walking. The LNK-Up Shoeless Option enables the foot to be lifted when not wearing shoes.

The band wraps around the midfoot, and fastens back on itself with hook/loop. The top part has a small loop to which the tensioning elastic is attached, in exactly same way as the pad, acting as a sort of replacement for the pad. Anatomically designed to fit the shape of the foot. Indications for use include;

- + Flaccid paralysis of the foot
- + Light to moderate spastic equinus foot

NOW  
AVAILABLE  
VIA **NHS**  
SUPPLY  
CHAIN

Allows for  
use of LNK-  
Up without  
shoes

LNK-Up  
Shoeless  
Option



Breathable  
micropore fabric  
ankle brace



Pad with metal  
hook to attach  
brace to shoe



Tensioning elastic  
with two hook/  
loop tips



Optional metal  
hooks to fit  
without shoeless  
option

Size	LNK-Up	NHSSC Code	Ankle Circ.	LNK-Up Shoeless Option - Left	NHSSC Code	LNK-Up Shoeless Option - Right	NHSSC Code	Foot Circ.
XS	LNK/02	GNA85751	12cm - 16cm	-	-	-	-	-
S	LNK/03	GNA85752	16cm - 20cm	LNKS/13*	GNA85748	LNKS/23*	GNA85755	17cm - 21cm
M	LNK/05	GNA85749	20cm - 24cm	LNKS/15	GNA85750	LNKS/25	GNA85747	21cm - 25cm
L	LNK/07	GNA85756	24cm - 28cm	LNKS/17	GNA85754	LNKS/27	GNA85753	25cm - 29cm

\*LNKS/13 and LNKS/23 are compatible with both LNK/02 and LNK/03 sizes.

Size	Replacement LNK-Up Elastic Attachment
XS	LNK/EA/02
S	LNK/EA/03
M	LNK/EA/05
L	LNK/EA/07

Size	Replacement LNK-Up Shoe Insert with Hook
Universal	LNK/SH
Size	Replacement LNK-Up Metal Hooks (pack of 4)
Universal	LNK/MH