


Paediatric Thumb Splint

Paediatric thumb abduction and opposition for patients with cerebral palsy, hemiplegia, indwelling thumb and other conditions such as congenital clasped thumb.

To size, measure the circumference of the distal thumb joint (A) and the distance on the palm from the distal thumb joint to the centre of the wrist crease (B). Refer to chart provided.



Size	Left	Right	A	B
1	MS-1L	MS-1R	3cm - 3.8cm	3cm - 4cm
2	MS-2L	MS-2R	3.4cm - 4cm	3.6cm - 5cm
3	MS-3L	MS-3R	3.7cm - 4.5cm	4cm - 6.3cm
4	MS-4L	MS-4R	4cm - 5cm	4.5cm - 6.9cm
5	MS-5L	MS-5R	4.3cm - 5.5cm	5.5cm - 7.9cm
6	MS-6L	MS-6R	4.8cm - 6cm	6.2cm - 9cm