

Mini Thumb Wrap

A one-piece low-profile thumb and wrist wrap designed to provide a gentle level of compression at the base of the thumb. This support assists with the treatment of conditions like osteoarthritis, fatigued joints and CMC and MCP instability or weakness.

Provides wrist and thumb support for an array of injuries, conditions, pre-and post-operations such as:

- + Arthritis
- + Carpel Tunnel
- + CMC and MCP Issues
- + Fractures
- + Ligament Injuries
- + Osteoarthritis
- + Thumb Fatigue



Size	Left	Right	Wrist Circ.
XS	MTW/12	MTW/22	11cm - 13cm
S	MTW/13	MTW/23	14cm - 15cm
M	MTW/15	MTW/25	16cm - 17cm
L	MTW/17	MTW/27	17cm - 18cm
XL	MTW/18	MTW/28	18cm - 19cm