

Comfort Cool™ Firm D-ring Wrist & Thumb



Multiple stays provide greater immobilisation of acute soft tissue injuries of the thumb and wrist. Made of 3.2mm perforated neoprene with a terry cloth liner that helps keep skin dry. Circumferential design with D-ring strap closures under double-layer stockinette makes this splint easy to put on.

Includes two dorsal stays and one volar wrist stay to prevent unwanted flexion and extension while allowing full use of the fingers. Thumb stay prevents unwanted motion of all three joints. Thumb tip can be trimmed for IP flexion. Volar stay is pre-bent to position the wrist in functional extension (approximately 15°).

- + Tendinitis
- + De Quervain's
- + Sprains & strains
- + Post surgical repair

Size	Left	Right	Wrist Circ.
XS	NC52992	NC52993	18cm - 20cm
S	NC52994	NC52995	20cm - 21cm
M	NC52996	NC52997	23cm >