

Comfort Cool[™] CMC Thumb Abduction Splint



Comfortable splint provides direct support for the thumb CMC joint while allowing full finger function. Microwave the thumb strap to soften the thermoplastic insert, then mould it to fit precisely along the thenar eminence crease. Thumb strap helps to abduct and stabilise the first metacarpal and open the thumb web space.

- + Thumb deformity
- + Contractures due to neurological conditions
- + Median nerve injuries
- + Arthritis (where the thumb is held abducted against the palm)

Size	Left	Right	DPC Circ.
S	NC79494	NC79495	15cm - 18cm
M	NC79496	NC79497	18cm - 20cm
L	NC79498	NC79499	20cm - 23cm