

Comfort Cool™ Wrist Wrap with Ulnar or Dorsal Support



The 1.6mm perforated neoprene supports without excessive bulk. Double-layer insert pocket protects and supports either the dorsal or the ulnar side of the wrist and helps prevent unwanted movement. Soft, terry cloth line feels good against the skin. Includes a 2.4mm thick thermoplastic insert offering firm support and restricting wrist movement depending on placement. Allows thumb and fingers to move freely. For dorsal support, simply order the splint size for the opposite hand. Anti-microbial and latex free.

- + Prevents ulnar deviation
- + Prevents wrist extension
- + Wrist instability
- + Ulnar drift problems
- + Triangular fibrocartilage

Size	Left - Regular	Right - Regular	Left - Long	Right - Long	Wrist Circ.
XS	NC91250	NC91251	NC91260	NC91261	15cm - 18cm
S	NC91252	NC91253	NC91262	NC91263	18cm - 20cm
M	NC91254	NC91255	NC91264	NC91265	20cm - 21cm
L	NC91256	NC91257	NC91266	NC91267	21cm - 23cm
XL	NC91258	NC91259	NC91268	NC91269	23cm >