

Comfort Cool™ Thumb Abductor Strap



Made of 1.6mm thick neoprene strap that provides a comfortable stretch for the finger or thumb MP's. Use to decrease thumb adductor tone or to position the thumb for optimal hand function. Achieve radial or palmar thumb abduction by relocating the thumb loop on the separate wrist strap. Can be used for trigger finger treatment with loop placed around MP and across dorsal hand to wrist strap. Fits either the left or right hand. Latex free. To size, measure circumference at the wrist crease. Small size is appropriate for some children.

- + Thumb spasticity
- + Tendinitis
- + Trigger finger

Size	Code	DPC Circ.
S	NC79520	11cm - 14cm
M	NC79521	14cm - 18cm
L	NC79522	18cm - 20cm