

Neoprene Thumb Spica

A Neoprene thumb spica featuring a removable flexible stay which restricts movement at the CMC and MCP joints, leaving the IP joint free.

Provides thumb CMC and MCP joint support for an array of injuries, conditions, pre- and post-operations such as:

- + Arthritis
- + Sprains and strains
- + Soft tissue injuries and conditions
- + Instability



Size	Left	Right	Wrist Circumference
S	NTS/13	NTS/23	14cm - 15cm
M	NTS/15	NTS/25	15cm - 16.5cm
L	NTS/17	NTS/27	16.5cm - 18cm
XL	NTS/18	NTS/28	18cm - 19cm