

Nuform Wrist Brace



An 18cm wrist brace made from ventilated elastic, which provides gentle compression and excellent breathability. This wrist brace features a vinyl pocket with an aluminium palmar bar for increased support and stabilisation, which can be shaped to provide a personalised fit. Available in beige or black.

This orthopaedic support offers a fantastic stabilisation at the wrist. By stabilising the joint, the wrist brace ensures that there is no excessive movement relative to the surrounding structures.

Provides wrist and thumb support for an array of injuries, conditions, pre-and post-operations such as:

- + Arthritis
- + Carpal tunnel syndrome
- + Inflamed/swollen tendons
- + Soft tissue injury
- + Degenerative irritation
- + Tendonitis
- + Worn cartilage
- + Wrist trauma



Size	Beige - Left	Beige - Right	Black - Left	Black - Right	Wrist Circ.
XS	NUF/12	NUF/22	NUF/12/BK	NUF/22/BK	11.5cm - 13cm
S	NUF/13	NUF/23	NUF/13/BK	NUF/23/BK	14cm - 15cm
M	NUF/15	NUF/25	NUF/15/BK	NUF/25/BK	16.5cm - 18cm
L	NUF/17	NUF/27	NUF/17/BK	NUF/27/BK	20cm - 21.5cm
XL	NUF/18	NUF/28	NUF/18/BK	NUF/28/BK	23cm - 24cm