

Oval-8 Finger Splint

Made of 1.6mm high-temperature polypropylene, Oval-8 splints can be adjusted with a heat gun for a small angle alteration. With a simple turn of the splint, each Oval-8 Finger Splint can be used to correct Swan Neck, Boutonniere and Mallet deformities, reduce Lateral Deviation and rest Trigger Finger.

Oval-8 splints feature a wide, flat band for comfortable pressure distribution and control. Each splint fits both a full and half ring size. Sizes correspond to jeweler's ring sizes, though correlation may vary due to swelling or joint changes.

The economical Oval-8 Kit includes 42 splints in sizes 2-14, a carrying case and instructions. Packaged with most splints in the more commonly used mid-range sizes. The Oval-8 Sizing Set includes one each of splint size 2-14 and the ring gauge. Splints sold in packs of 5.

Size	Beige	Clear
Complete Kit	P1008-K	-
Sizing Set	P1008-S	P1108-S
2	P1008-5-02	P1108-5-02
3	P1008-5-03	P1108-5-03
4	P1008-5-04	P1108-5-04
5	P1008-5-05	P1108-5-05
6	P1008-5-06	P1108-5-06
7	P1008-5-07	P1108-5-07
8	P1008-5-08	P1108-5-08
9	P1008-5-09	P1108-5-09
10	P1008-5-10	P1108-5-10
11	P1008-5-11	P1108-5-11
12	P1008-5-12	P1108-5-12
13	P1008-5-13	P1108-5-13
14	P1008-5-14	P1108-5-14
15	P1008-5-15	P1108-5-15



Complete Kit