

Procool D-Ring Wrist Splint

Low profile support with circumferential slip-on design and D-Ring fastening to allow easier one handed application. Aluminium palmar bar can be removed for further adjustment. Dorsal stockinette lining facilitates application and cushions strap closures. Design of the Standard and Firm braces allow full finger and thumb function whilst limiting wrist flexion.

The Firm splint has the addition of adjustable dorsal stays to reduce unwanted flexion and extension and improve overall immobilisation of acute soft tissue injuries. This brace can be trimmed if required.

Provides wrist and thumb support for an array of injuries, conditions, pre-and post-operations such as:

- + Arthritis
- + Sprains & strains
- + Carpal Tunnel Syndrome



Size	Code - Left	Code - Right	Code - Left (Firm)	Code - Right (Firm)	Wrist Circumference
XS	PDR/12	PDR/22	PDRF/12	PDRF/22	< 15cm
S	PDR/13	PDR/23	PDRF/13	PDRF/23	15cm - 17cm
M	PDR/15	PDR/25	PDRF/15	PDRF/25	17cm - 19cm
L	PDR/17	PDR/27	PDRF/17	PDRF/27	19cm - 22cm
XL	PDR/18	PDR/28	PDRF/18	PDRF/28	22cm >