

Paediatric Thumb Spica

A thumb orthosis which provides adjustable support to the thumb carpometacarpal (CMC) and metacarpal phalangeal (MCP) joint. A flexible removable stay restricts movement at the CMC and MCP joints, leaving the interphalangeal (IP) joint free.

Provides wrist and thumb support for an array of injuries, conditions, pre-and post-operations such as:

- + Sprains and strains
- + Soft tissue injuries and conditions
- + Instability
- + Arthritis



Size	Left	Right	Wrist Circ.
XS	PLTS/12	PLTS/22	12cm - 13cm
S	PLTS/13	PLTS/23	13cm - 14cm
M	PLTS/15	PLTS/25	14cm - 15cm
L	PLTS/17	PLTS/27	15cm - 17cm
XL	PLTS/18	PLTS/28	17cm - 18cm