

# Paediatric Wrist / Thumb Brace

Supporting and protecting the wrist after an injury occurs is crucial to the healing and rehabilitation process, especially in children who are active, often unaware and on-the-go more than adults. The Paediatric Wrist/Thumb Splint grants your child an optimal balance between comfort and protection while they recover from sprains, strains, instability fractures and a range of other soft-tissue injuries.

Provides wrist and thumb support for an array of injuries, conditions, pre-and post-operations such as:

- + Arthritis
- + Sprains and strains
- + Soft tissue injuries and conditions
- + Instability
- + Post fracture management



Lycrafleece material is designed to be tough, durable and comfortable, while remaining breathable to help with temperature regulation.

Size	Left	Right	Wrist Circ.
XXS	PLWTB/11	PLWTB/21	11cm - 12cm
XS	PLWTB/12	PLWTB/22	12cm - 13cm
S	PLWTB/13	PLWTB/23	13cm - 14cm
M	PLWTB/15	PLWTB/25	14cm - 15cm
L	PLWTB/17	PLWTB/27	15cm - 17cm
XL	PLWTB/18	PLWTB/28	17cm - 18cm