

Neorpene Wrist / Thumb Brace

A wrist thumb orthosis with malleable aluminium thumb stay to hold the thumb firmly in abduction, supporting the thumb MP joint and restricting movement at the CMC joint.

Provides wrist and thumb support for an array of injuries, conditions, pre-and post-operations such as:

- + Arthritis
- + Sprains and strains
- + Soft tissue injuries and conditions
- + Instability
- + Post fracture management



Size	Black - Left	Black - Right	Wrist Circ.
S	PNWT/13/BK	PNWT/23/BK	11.5cm
M	PNWT/15/BK	PNWT/25/BK	14.5cm
L	PNWT/17/BK	PNWT/27/BK	16cm
Size	Royal Blue - Left	Royal Blue - Right	Wrist Circ.
S	PNWT/13/RB	PNWT/23/RB	11.5cm
M	PNWT/15/RB	PNWT/25/RB	14.5cm
L	PNWT/17/RB	PNWT/27/RB	16cm
Size	Red - Left	Red - Right	Wrist Circ.
S	PNWT/13/RD	PNWT/23/RD	11.5cm
M	PNWT/15/RD	PNWT/25/RD	14.5cm
L	PNWT/17/RD	PNWT/27/RD	16cm