

Pro-Rheuma Thumb Spica

Made from a breathable high compression material with microfleece lining and metal stay for added rigidity. The Pro-rheuma Thumb Spica features a removable semi-flexible thumb stay that helps stabilise CMC and MCP joints.

Provides wrist and thumb support for an array of injuries, conditions, pre-and post-operations such as:

- + Arthritis
- + Sprains and strains
- + Soft tissue injuries and conditions
- + Instability



Size	Left	Right	Wrist Circ.
XS	PRHS/12	PRHS/22	13cm - 14cm
S	PRHS/13	PRHS/23	14cm - 16cm
M	PRHS/15	PRHS/25	16cm - 18cm
L	PRHS/17	PRHS/27	18cm - 20cm
XL	PRHS/18	PRHS/28	20cm - 22cm