

Promedics Shoulder Immobiliser

A shoulder orthosis designed to provide immobilisation. The sling supports the arm from wrist to humerus, with Y-tab adjustable straps for the shoulder and body, holding the arm close to the torso. The thumb loop limits the rotation of the shoulder, keeping it secure.

The universal design fits both the right and left arm. Suitable for a range of injuries/conditions where restriction of shoulder and/or elbow movement is required, such as:

- + Dislocation
- + Strains and sprains
- + Pre and post operatively
- + Rehab stages
- + Brachial plexus or peripheral nerve



Size	Code	Forearm Length
XS	PSI/02	< 33cm
S	PSI/03	33cm - 37cm
M	PSI/05	37cm - 40cm
L	PSI/07	40cm - 43cm
XL	PSI/08	43cm >