

Procool Wrist Thumb Support

An adjustable wrist thumb wrap offering good support around the wrist, MCP and CMC joints. Adjustable thumb closure with excellent shaping around the palm.

Procool is a range of lightweight and breathable splints producing excellent compression and support for a wide range of conditions and injuries. Made with tiny perforations in the Neoprene for comfortable wear and a terry cloth lining that feels good against the skin. Procool offers you just the fit you want - excellent stability without excessive bulk.

Provides wrist support for an array of injuries, conditions, pre-and post-operations such as:

- + Arthritis
- + Sprains and strains
- + Soft tissue injuries and conditions
- + Instability



Size	Left	Right	Wrist Circumference
XS	PWTS/12	PWTS/22	12cm - 14cm
S	PWTS/13	PWTS/23	14cm - 16cm
M	PWTS/15	PWTS/25	16cm - 18cm
L	PWTS/17	PWTS/27	18cm - 20cm
XL	PWTS/18	PWTS/28	20cm - 22cm