

# Orfilight<sup>®</sup>

## Dorsal Wrist Cock-Up Splint

An excellent alternative to volar wrist cock-up orthosis as this design allows for increased function while maintaining the wrist in extension.

The orthosis design is suitable as a base or the attachment of low profile outriggers.



Code	Style	Size
35818/L	Beige • 3.2mm Non perfo	XS
35820/L	Beige • 3.2mm Non perfo	S
35821/L	Beige • 3.2mm Non perfo	M
35822/L	Beige • 3.2mm Non perfo	L
35820Z/L/NS	Black • 3.2mm Non perfo	S
35821Z/L/NS	Black • 3.2mm Non perfo	M
35822Z/L/NS	Black • 3.2mm Non perfo	L