

# Air X Wrist Thumb

A two-piece wrist thumb brace with longitudinal seam conforms perfectly to the wrist, hand and forearm. Aluminium palmar bar and thumb stay provides excellent support to wrist and first CMC & MCP joints. The thumb stay can be easily adjusted for the required degree of abduction and the thumb closure is easily adjusted for a comfortable fit

Provides wrist and thumb support for an array of injuries, conditions, pre-and post-operations such as:

- + Arthritis
- + Sprains and strains
- + Soft tissue injuries and conditions
- + Instability





The unique Air X fabric transports moisture away from the body. The layered fibre structure promotes a moisture wicking process, resulting in a dry comfortable feeling so the skin always feels fresh.



Can be easily adjusted for a personalised fit unique to the user, with a longitudinal seam that conforms perfectly to the wrist.



Aluminium palmar bar and thumb stay provide excellent joint support, palmar bar can be removed if necessary or desired.

Code (Left)	Code (Right)	Size
AWT/12	AWT/22	XS
AWT/13	AWT/23	S
AWT/15	AWT/25	M
AWT/17	AWT/27	L
AWT/18	AWT/28	XL